

CLINISAN for Commercial Cleaning

CLINISAN is a simple, accurate and safe solution. It's versatility allows for streamlining of products for any commercial cleaning operation.



Benefits

- **Less risk** to your staff when diluting disinfectant
- **Less storage room** required
- **pH of 5.5-7** - Close to skin pH and kinder on staff and patients
- **Less plastic waste** creation and therefore less waste processing costs
- **Lower risk Dangerous Goods Class** than typical Chlorine compounds
- **Less corrosive** than acidic, alkali and other chlorine products
- **Accurate dilution** - over or under dilution risk reduced

Dilution	Concentration (ppm)	Tablets Dilution (1g tablet per)	Tablets Dilution (5g tablet per)	Use
Hospital Grade Disinfection	1000	300 ml	1.5 L	Spray / Wipe / Mop
Critical Toilet and Bathroom	1000	300 ml	1.5 L	Spray / Wipe / Mop
Day-to-Day Use	150	2 L	10L	Spray / Wipe / Mop
Alternative Day-to-Day Use	400	1 tablet per spray bottle		Spray / Wipe
No-Rinse Food Surface (MPI's C42)	150	2 L	10 L	Spray / Wipe
Toilet Cleaner		1 tablet per bowl		Once dissolved, use as normal
Mops and Dishcloth Soak	60	5 L	25 L	Soak for 1 hour



Bactericidal
99.999%
5 min contact



Norovirus
99.99%
5 min contact



Results at 1000 ppm

Tested to TGA's Hospital Grade Standard



Fungicidal
99.99%
5 min contact



Coronavirus
99.999%
5 min contact

Health and Safety

CLINISAN's active ingredient was tested in a confined space compared to bleach. The airborne chlorine was measured in this assessment to see the exposure to users of each disinfectant

	Concentration (ppm)	ACGIH Standard	Peak Detected Airborne Chlorine (ppm)	% of samples over Lower Detection Limit
Bleach	0.2%	0.4 ppm	0.5	31%
	0.5%		0.7	65%
NaDCC (Clinisan)	0.2%		0.1	2%
	0.5%		0.1	3%

Lower Detection Limit of 0.1 ppm

ACGIH - American Conference of Governmental Industrial Hygienists - **The strictest standard to be found**

NaDCC is was found to have between 2% and 5% measurements over the Lower Detection Limit compared to Bleach's 31% and 65%, also a much lower peak of 0.1 ppm compared to 0.7 ppm. Showing a much reduced exposure to airborne chlorine for your staff and children

Biodegradability and Sustainability

NaDCC has been found to be biodegradable under real world conditions producing beneficial environmental products.

Clinisan is stored in a fully recyclable tub with no individual packaging.